



GIBSON LABS TRI-VALENT™ SWAB

THREE CONTROLS • ONE SWAB

Gibson Labs offers positive and negative control swabs for use with the BD Affirm™ VPIII Microbial Identification System DNA Probe Detection Kit. The Tri-Valent™ Swab functions as a tri-analyte swab that can be processed following the exact methods for a patient specimen.

Each swab includes whole stabilized cells of *Trichomonas vaginalis*, *Gardnerella vaginalis*, and *Candida albicans* and is easy to use with extraction vials.

Cat. # TVS-01 Tri-Valent™ Swab (Positive Control)

Positive control for BD's Affirm VPIII Microbial Identification for *Trichomonas*, *Gardnerella*, and *Candida*

Cat. # TVS-02 Tri-Valent™ Swab (Negative Control)


Negative control for BD's Affirm VPIII Microbial Identification for *Trichomonas*, *Gardnerella*, and *Candida*

Cat. # VV-01 Tri-Valent™ Swab (Validation Set)

Randomly mixed set of 22 positive and negative controls for BD's Affirm VPIII Microbial Identification for *Trichomonas*, *Gardnerella*, and *Candida*

Cat. # CV-01 Tri-Valent™ Swab (Calibration Panel)

Set of 6 swabs containing high and low values for BD's Affirm VPIII Microbial Identification for *Trichomonas*, *Gardnerella*, and *Candida*



Visit www.gibsonlabs.com or call 800.477.4763 to place your order!

1040 Manchester Street Lexington, KY 40508 • Phone: 859-254-9500
800-477-4763 • Fax: 859-253-1476
www.gibsonlabs.com